


Application of CoolFit on the Inner Thigh



Step 1. Assess the treatment area while the patient is standing with the untreated leg slightly raised.



Step 2. Mark the peak of the bulge and treatment site, using the CoolFit treatment template as needed.



Step 3. Place the gelpad over the treatment area, ensuring the peak of the bulge is covered.



Step 4. With one hand supporting the leg, place the CoolFit applicator over the treatment area.



Step 5. During the treatment, the patient can be seated with the leg slightly bent. A pillow can be used to support the lower leg.



Step 6. After treatment, massage the treatment area for 2 minutes